

TECHNICAL DATA SHEET
O2 FRESH™ Nasal Cannula Cover

SANITIZABLE AND REUSABLE

Intended Use: Hygienic storage of nasal cannula prongs; patient comfort and convenience, hang tubing out of way and declutter bedside.

Part description: Light gray with semi-matt textured finish.

Not an FDA classified device.



Description	SKU	UOM
Nasal Cannula Cover (1)	THG202101	Each
Nasal Cannula Cover (case of 50)	THG202101	Case

Features	Universally attachable handle, ventilation holes, opening with secure, non-damaging internal clamp.	Material Certifications	Prop 65 compliant, UL 94 certification for flammability, REACH RoHS
Areas of Use	Hospitals, senior living facilities, skilled nursing facilities, post-acute care, homecare	DEHP	No
Sterile	No	Biocompatibility	USP Class I-VI and the U.S. FDA guidance for Intact Skin Surfaces
Patient Population	Adult	Packaging Options	Individually packaged 50/case

Cleaning Techniques

- Top shelf dishwasher safe.
- Clean with mild detergent.

Sterilization Techniques

- Boil in water for 5 minutes. **DO NOT MICROWAVE.**
- Run through dishwasher, setting the water temperature to 158° or more for 30 minutes. Top shelf only.
- Soak it in 70% isopropyl alcohol for 5 minutes, rinse using sterile water, and air dry before use.
- Soak it in 3% hydrogen peroxide for 30 minutes, rinse using sterile water and air dry before use.

See website for list of compatible adult low flow and high flow nasal cannulas.

Warnings: Not a toy; keep out of reach of children. Not microwave safe.